
### Snoring Detection Based on A Stretchable Strain Sensor


Zhengwei Li<sup>1, 2</sup>, Qingkun Song<sup>1</sup>, Long Cheng<sup>2, 3</sup>, Min Tan<sup>2, 3</sup>

 School of Automation, Harbin University of Science and Technology, Harbin 150080, China.
State Key Laboratory of Management and Control for Complex Systems, Institute of Automation, Chinese Academy of Sciences, Beijing 100190, China.
University of Chinese Academy of Sciences, Beijing 100049, China.
















# Snoring detection







# Introduction- risks of snoring

Influence on daily lives

- sleep deprivation;
- daytime drowsiness;
- irritability;
- □ lack of focus;
- □ decreased libido;

- Risks of suffering from diseases
- 67% greater chance of suffering from

stroke

**3**4% increased risk of having heart


attack



Source: [1] http://news.bbc.co.uk/2/hi/health/7272651.stm Source: [2] https://goodmorningsnoresolution.com/pages/how-it-works Obstructed Airways Cause Most Snoring

# Introduction- detection methods

#### Polysomnography (PSG)




- Characteristics of polysomnography
- A standard method of monitoring the sleeping quality in the literature;
- □ Acquire various signals (EEG, EMG, ECG, EOG, Respiration condition);
- □ Require the experts among analysis of the measurement result;
- **D** Expensive, not suitable for snores to use at home.

[2]

# Introduction- detection methods

#### Audio sensor





- Characteristics of using audio data for snoring detection
- □ Low price and small volume;
- □ Prone to be affected by the environment noise;

the machine-learning-based detection method requires a large amount of computation resources.

# **Introduction- detection methods**

Piezoelectric vibration sensor

□ Signal strength depends on whether the piezoelectric sensor is tightly

attached to the snorer's skin.

- Nasal cannula -----

bring discomfort to the snorers;

easily affected by nasal blockage, mouth breathing and the movement of the

cannula in the nares.

Is there a natural and friendly way to detect snoring and avoid issues mentioned before?










# **B** Shoring detection




# 2 Strain sensor: assemble processing



Fabrication process of stretchable strain sensor.





The sensor is low-cost, light and portable, and can be use at home.



Light intensity received by silicon photocells

$$E_l = Echv \frac{\Delta L}{L_0}$$

Where E is the light intensity emitted by the LED, and v is the Poisson's ratio of the rubber;


Sensor output voltage


$$\Delta U = K_d K_i R_c E ch v_{\overline{L_0}}^{\Delta L}$$

Where  $K_d$  is the amplification factor of the amplifier,  $R_c$  is the series resistance of the circuit, and  $K_i$  is the slope of the current generated by the silicon photocell and the intensity of the received light.


#### Sensor sensitivity

$$GF = \frac{\Delta U}{\Delta L/L_0} = K_d K_i R_c Echv$$





**2** Strain sensor: performance test



Strain sensor has a high accuracy and a good linearity.





## Snoring detection: algorithm

#### Snoring detection principle

Amplitude-based detection principle:
(1) Determine the amplitude threshold A<sub>light</sub> to distinguish in normal sleeping or not ;
(2) Determine the fluctuation lasting time T<sub>last</sub> to distinguish in Short-term actions or continuous actions;

(3) Determine the amplitude threshold

 $A_{\text{snore}}$  to distinguish in snoring or deep breathing.

Slope-based detection principle: Determine the slope  $P_{\text{max}}$  and  $P_{\text{min}}$  to distinguish in snoring or not ; Algorithm 1 Sleeping action recognize

- 1: while (1) do
- 2:  $t = 0, a = 0, p = 0; \bigcirc$
- 3: obtain 10 cycles of fluctuating data D;
- 4: unbiased processing of fluctuating data D E(D);
- 5: search the local maximum values  $D_m$  and local minimum values  $D_n$ ;
- 6: while i < 10 do
- 7: calculate the time difference between two adjacent extremum points  $T(i) = t(D_m(i)) - t(D_m(i-1));$

calculate the relative amplitude of each cycle  $A(i) = D_m(i) - D_n(i);$ 

calculate the slope of signal rising edge  $P(i) = \frac{D_m(i) - D_n(i)}{t(D_m(i)) - t(D_n(i))}$ .

- 10: if  $A(i) > A_{\text{light}}$  then
  - if  $A(i) > A_{\text{snore}}$  then
  - a = a + 1;end if
  - if  $A(i) A(i-1) < 0.5 \times A(i)$  then
- 15: t = t + 1;
- 16: else

11:


12:

13:

14:

- 17: output t and t = 0;
- 18: end if
- 19: if  $P(i) > P_{\min}$  and  $P(i) < P_{\max}$  then
- 20: p = p + 1;
- 21: end if
- 22: end if
- 23: end while
- 24: if  $a > A_{\text{snore}}$  and  $\max(t) > T_{\text{last}}$  and p > 5 then
- 25: snoring happened during sleeping.
- 26: end if
- 27: end while


# **Snoring detection:** experiments



Five subjects are employed to test the proposed snoring detection method.

# Snoring detection: result


Audio data and strain sensor data are recorded when snoring occurs during sleeping.



It can be observed that strain sensor is sensitive to snoring.

# Snoring detection: result

Audio data and strain sensor data is recorded when snoring and swallowing occur during sleeping.




# Snoring detection: conclusion

This MOOP proposes a new method to detect snoring based on a selfdeveloped stretchable strain sensor.

#### Conclusions

- This proposed detection method can avoid the influence of environment noise during the detecting process;
- This sensor only applies light stress to the snorer's skin and obtains a satisfactory detection result;
- This detection device is small, light and cheaper, it can be used by snorers at home;
- The algorithm of detecting snoring/non-snoring actions was proposed and verified by experiments.









Develop a more convenient way to wear strain sensor;
Machine learning methods can be used to detect the snoring, swallowing and turning over actions and therefore monitor the sleeping quality;
Other interesting applications based on the self-developed strain sensor are to be investigated.

# Thanks for your attention